



Rosen Knee Revision Tibial Retractor

Designed modification by Adam Rosen, DO of original design by Christopher M. Meckel, MD

Designed for total knee and revision total knee replacements using posterior stabilized knee components, the posterior build-up on the retractor allows the surgeon to more easily translate the tibia forward for better visualization after femoral notch preparation



PRODUCT NO:
2830
Overall Length: 12.5" (31,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362