

Rotating Offset Handle Hex Driver

*Offset shaft and smooth spin handle allow for
a rapid crank action when desired*

Large hex driver for 6.5 mm and 4.5 mm diameter screws.
Especially helpful in insertion and removal of long screws.

PRODUCT NO:

7241

Overall Length: 13.5" (34,3 cm)

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362