

# Sarraf Toothed Curettes

Designed by Khaled Sarraf, MD

*Forward, straight, and reverse bent toothed curettes designed to aid in aid in all types of joint arthroplasty surgery, especially in scraping any articular chondral islands within the acetabulum during THA preparation*

- ▶ Can also be used for the femoral canal in cemented and uncemented THA
- ▶ Valuable aid in revision arthroplasty (hip, knee, shoulder and ankle) for cement curettage
- ▶ Useful tool in hip and knee primary arthroplasty as well as shoulder, elbow and ankle arthroplasty procedures



PRODUCT NO'S:
5174-00 [Set]
<b>Set Includes / Available Separately:</b>
5174-F [Forward Toothed Curette] Overall Length: 11.5" (29,2 cm) Handle Length: 5.5" (14 cm) Angled Down: 30°
5174-R [Reverse Toothed Curette] Overall Length: 11.5" (29,2 cm) Handle Length: 5.5" (14 cm) Angled Up: 30°
5174-S [Straight Toothed Curette] Overall Length: 11.5" (29,2 cm) Handle Length: 5.5" (14 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362