


# Screw/Pin Removal Locking Pliers

Unique jaw designed to solidly grip and clamp onto a screw head, broken screw, or pin for removal



**NEW Small**  
Jaw End & Bite

Designed to securely grab pins as small as 1.4 mm (.055") up to 2.4 mm (.095")

**Standard**  
Jaw End & Bite

Designed to securely grab larger pins, screw heads, or broken screws




Standard

**New!**

Small

Reduced jaw size for smaller screws, pins and incisions

PRODUCT NO'S:	
S0142 [Standard]	Overall Length: 7.875" (20 cm) Jaw Width at End: 4 mm
S0142-01 [Small]	Overall Length: 7.875" (20 cm) Jaw Width at End: 4 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** <sup>®</sup>



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.