



# Seymour ACL Graft Advancer

Designed by Scott Seymour, MD

*Designed to facilitate the passage and tensioning of an ACL graft into the femoral and tibial tunnels*

A loop is tied in the prepared graft's passing sutures and the device is used to pull the graft into the tunnels, then to tension the fixation.



**PRODUCT NO:**

1117

Overall Length: 4.35" (11,1 cm)

Handle Width: 4" (10,2 cm)

Hook Width: 19,5 mm Outside, 13,5 mm Inside

Hook Depth: 25 mm

Hook Diameter: 3 mm



USA MADE

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**