

Sherman Patella Tendon Harvest Retractor



Designed by Mark Sherman, MD

Designed to help improve exposure and lessen the incision necessary to harvest a patella tendon graft during anterior cruciate ligament bone–patella tendon–bone (BTB) reconstruction

PRODUCT NO:

4691

Overall Length: 8.5" (21,6 cm)
Depth from Bend: 3.25" (8,3 cm)
Blade Width: .875" (22 mm)
Prong Length: 10 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED