



Shouldered Pin Retractors

Designed by M. Jake Hamer, MD

Useful in fracture and reconstructive cases

Can be used for opening wedge osteotomies, obtaining length or distraction to correct a foreshortened malunion, or in trauma cases (calcaneus fractures, tibial plateau fractures, pilon fractures). With three different length shouldered pins, the retractor can be set at different depths to provide the necessary retraction and/or distraction.

Can be used with any compressor and/or distractor that can accept 1/8" (3,2 mm) pins.



PRODUCT NO'S:
Sold In Pairs
4694-01 [Long] Overall Length: 67 mm End-to-shoulder: 22 mm Pad: 13 mm x 15 mm
4694-02 [Medium] Overall Length: 67 mm End-to-shoulder: 31 mm Pad: 13 mm x 15 mm
4694-03 [Short] Overall Length: 67 mm End-to-shoulder: 40 mm Pad: 13 mm x 15 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED