

Sidhu Tibia Clamp

Designed to be used to securely grasp and remove an entire tibial wedge

Designed by
Kuldeep Sidhu, MD

Product No: 3643
Overall Length: 10.25"
Pads: 60mm x 30mm
Spike Length: 7.5mm



The tapered lower pad slides under the cut tibial wedge without first having to use wedges, then, clamping allows the spikes in the upper pad to securely grasp the entire tibial wedge for easy removal.



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

