

Sierra OrthoLucent™ EVA Pelvic Osteotomy Retractor

Designed by Rafael J. Sierra, MD

Designed to help with retraction of the inner pelvis for direct visualization of the inner pelvis prior to iliac osteotomy, the retractor is made of a lightweight carbon fiber PEI composite material



The OrthoLucent™ carbon fiber PEI material is radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

PRODUCT NO:

4541

Overall Length: 10.7" (27.2 cm)

Blade Width at Widest: 24,5 mm

Blade Width at End: 10 mm

MADE EXCLUSIVELY
FOR INNOMED IN
SWITZERLAND



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED