


Sierra OrthoLucent™ Soft Tissue Retractor

Designed by Rafael J. Sierra, MD

Radiolucent retractor designed for soft tissue protection of lateral muscles during pelvic osteotomy surgery

Manufactured of delrin and aluminum.

PRODUCT NO:
4849
Overall Length: 10.25" (26 cm)
Handle Length: 9" (22.9 cm)
Blade Depth: 4" (10.2 cm)
Blade Width: .75" (19 mm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED [®]

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.