




## Slavitt Phalangeal Forceps

Designed by Jerome Slavitt, DPM

*Designed to enable the surgeon to provide joint distraction and stability during prosthesis placement at the base of the proximal phalanx of the lesser digits*



Helps to distract the joint and hold the bone, allowing easier access to the base. Can also be used for digital fusions to hold bones better for drilling and cutting applications.

**PRODUCT NO:**

1163

Overall Length: 6"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362