

# Stulberg Sliding Bolster

Designed by S. David Stulberg, MD

*Helps eliminate the need for a sand bag during total knee surgery*


The base plate is attached to the table and the sterile sliding bolster is placed on top of the sterile drape. The bolster can be adjusted for different angles of knee flexion during surgery.

**PRODUCT NO'S:**

**2730**

Base Dimensions:

20" x 10.5" (50,8 cm x 26,7 cm)



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**