



Soft Impact Mallets with Easy Grip Handles

Provides shock-absorbing force

Designed to have a shock-absorbing force, providing less bounce or wasted force. The Mallet is filled with a shock-absorbing media and has a flat striking surface to keep the Mallet centered on an instrument.

The comfortable Easy Grip Handle is made of a textured silicone that helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip. The bottom can also be used to tap an implant in place.

PRODUCT NO'S:

7820 [2 lbs. (.907 kg) Standard]
Overall Length: 10.5" (26.7cm)
Handle Length: 5" (12.7cm)
Head Width: 3.5" (8.9cm)
Head Diameter: 1.375" (3.5cm)

7832 [2 lbs. (.907 kg) with Delrin End]
Overall Length: 10.5" (26.7cm)
Handle Length: 5" (12.7cm)
Head Width: 3.5" (8.9cm)
Head Diameter: 1.375" (3.5cm)

7837 [3 lbs. (1.36 kg) Standard]
Overall Length: 11" (27.9cm)
Handle Length: 5" (12.7cm)
Head Width: 3.5" (8.9cm)
Head Diameter: 1.875" (4.8cm)



Soft Impact Mallet with Weidman Silicone Handle

Handle designed by Kevin Weidman, MD

PRODUCT NO:
7821 [2 lbs. (.907 kg)]
Overall Length: 10.625" (27cm)
Grip Length: 5.5" (14cm)
Head Width: 3.5" (8.9cm)
Head Diameter: 1.375" (3.5cm)



Scan to
Launch Our
Website



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362