



### Easy Grip Textured Soft Silicone Handle



Comfortable grip helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip.

#### PRODUCT NO'S:

7820 [2 lbs. Standard]

Weight: 2 lbs. (.907 kg)

Overall Length: 10.5" (26.7 cm)

Handle Length: 5" (12.7 cm)

Head Width: 3.5" (8.9 cm)

Head Diameter: 1.375" (3.5 cm)

7821 [2 lbs. With Weidman Handle]

Weight: 2 lbs. (.907 kg)

Overall Length: 10.625" (27 cm)

Grip Length: 5.5" (14 cm)

Head Width: 3.5" (8.9 cm)

Head Diameter: 1.375" (3.5 cm)

7832 [2 lbs. With Delrin End]

Weight: 2 lbs. (.907 kg)

Overall Length: 10.5" (26.7 cm)

Handle Length: 5" (12.7 cm)

Head Width: 3.5" (8.9 cm)

Head Diameter: 1.375" (3.5 cm)

7837 [3 lbs. Standard]

Weight: 3 lbs. (1.35 kg)

Overall Length: 11" (27.9 cm)

Handle Length: 5" (12.7 cm)

Head Width: 3.5" (8.9 cm)

Head Diameter: 1.875" (4.8 cm)

#### Delrin Head Replacements for 7832:

7832-HEAD01 [.5" Stud] Single

7832-HEAD02 [.5" Stud] 3-Pack

7832-HEAD03 [.875" Stud] Single

7832-HEAD04 [.875" Stud] 3-Pack

Weidman handle designed by Kevin Weidman, MD



Replacement Delrin Heads

# Soft Impact Mallets with Easy Grip Handles

## Provides shock-absorbing force

Designed to have a shock-absorbing force, providing less bounce or wasted force. The mallet is filled with a shock-absorbing media and has a flat striking surface to keep the mallet centered on an instrument.

The comfortable Easy Grip handle is made of a textured silicone that helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip. The bottom can also be used to tap an implant in place.

The mallet with delrin head features a replaceable delrin head.



Replaceable  
Delrin Head

Weidman  
Silicone Handle

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

USA MADE

103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362