

Star Bit Driver Set

Set consists of four star bits – T10, T15, T20, & T25, a handle which accommodates any of the above bits, and a sterilization case. The drive end (A/O) is designed for easy and quick engagement with the universal instrument handle. The ergonomic, modular handle has two connection points, allowing for both straight and T-handle orientations.



PRODUCT NO'S:	
5194-00	[4 Star Bits w/Handle & Case]
5194-01	[4 Star Bits w/Case only]
Also sold individually:	
S0113	[Universal 4" Handle]
5194-10	[T10 with A/O End]
5194-15	[T15 with A/O End]
5194-20	[T20 with A/O End]
5194-25	[T25 with A/O End]
9003	[Case]

**PROUDLY
MADE
IN THE
USA**



Set in Storage Case

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362