


# Stowell Modified Posterior Acetabular Retractor

Designed by R.L. Stowell, MD

*Designed to be placed along the posterior rim of the acetabulum to facilitate exposure and acetabular preparation*



**PRODUCT NO:**  
7330  
Overall Length: 14.5" (36,9 cm)  
Handle Length: 5.25" (13,3 cm)  
Blade Depth: 4" (10,4 cm)  
Blade Width: 2.75" (7 cm)  
Prongs: 8 mm Wide x 6 mm Long



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

1.800.548.2362