



Strayer Retractor

Designed by Irvin Oh, MD

A lamina spreader with long thin blades designed to retract the soleus muscle and soft tissue for isolation and exposure of the gastrocnemius fascia for release



PRODUCT NO:
1869
Overall Length: 9.25" (23,5 cm)
Blade Length: 3.5" (8,9 cm)
Blade Width: .6" (1,5 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362