



Stulberg Incision Close Gelpi & Blade Set

Designed by S. David Stulberg, MD

Designed to help expose difficult to visualize areas at the end of incisions



PRODUCT NO'S:	
4269-00	[Set – 1 Gelpi & 1 Blade]
Also available Individually:	
4269-01	[Gelpi]
Overall Length: 7.25" (18,4 cm)	
Maximum Spread Width: 3.5" (8,9 cm)	
4269-02	[Blade]
Overall Length: 5.5" (14 cm)	
Blade Width: 1" (2,54 cm)	
Blade Bend-Back Angle: 130°	



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362