


Subscapularis Spreader

Reaches deep to help split the subscapularis in a Jobe approach

Also used for retracting a split deltoid in mini rotator cuff repairs.



PRODUCT NO:
T1005 [Standard] Overall Length: 7" (17,8 cm)
T1005-L* [Long] Overall Length: 9.25" (23,5 cm)

MADE FOR INNOMED IN
GERMANY

* PROUDLY
MADE IN
THE USA

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362