


Superior Retractor

Used for retraction around the acetabulum, can be self retaining with the use of 1/8" (3.2 mm) pins



PRODUCT NO:
6027
Overall Length: 12.875 (32,7 cm)
Depth from Bend: 5.25" (13,3 cm)
Blade Width: 25,4 mm
Holes for Pins Up to: 1/8" (3,2 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362