

# Sure Grip Soft Tissue Grasper

Designed by  
Andrew Glassman, MD

**Enables the surgeon to securely grasp soft tissue structures within the knee**

Incorporates a spike into its upper jaw with a matching recess in the lower jaw, enabling the surgeon to securely grasp soft tissue structures within the knee. Particularly useful for grasping the posterior horn of either the medial or lateral meniscus. Also useful when excising the cruciate ligaments, capturing loose bodies, holding the retinaculum during patellar preparation, and grasping the capsule during wound culture.



Product No's:  
3645-01 [5"]  
3645-02 [7"]  
3646-02 [7" with Lock] **Shown NEW!**  
3645-03 [9"] *Shown*

NEW STYLE  
AVAILABLE

**Now available  
with a lock!**  
(7" size only)

**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

**Phone** 800 548 2362  
912 236 0000  
**Fax** 912 236 7766  
**Internet** www.innomed.net  
**eMail** info@innomed.net

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

