

Swanson Elevator


Designed by Richard Ferkel, MD

Angular design helps to go around bone for retraction and elevation – especially useful in small bone surgery of the hand/wrist and foot/ankle

PRODUCT NO:

1644

Overall Length: 6.375" (16,2 cm)
Blade Depth: .75" (1,9 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362