

# T-Handle Femoral Canal Finders


*Designed to sound the femoral canal prior to stem broaching, especially useful to help start the broach path during the direct anterior approach*

## Rockowitz T-Handle Femoral Canal Finder Rasp

Designed by Neal L. Rockowitz, MD

**PRODUCT NO:**

4990 [Rockowitz Rasp]  
Overall Length: 9" (22,9 cm)




**ORIGINAL  
DR. ROCKOWITZ DESIGN  
Topside Rasp**

Rasp on curve topside and sides,  
smooth on underside

## Modified T-Handle Femoral Canal Finder Rasp

**PRODUCT NO:**

4989 [Modified Rasp]  
Overall Length: 9" (22,9 cm)



**MODIFIED DESIGN  
Underside Rasp**

Rasp on curve underside and sides,  
smooth on topside

## T-Handle Femoral Canal Finder - Smooth

**PRODUCT NO:**

4990-03 [Smooth]  
Overall Length: 9.385" (24,4 cm)



**SMOOTH DESIGN**

No rasp – smooth underside,  
sides, and topside

ISO 13485:2016

Scan to  
Launch Our  
Website



**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**