




TKA Gap Assessment Gauge Assembly

Designed by Michael Radon

Universal design allows the gauge to be used without the removal of trials to help determine if a 1 or 2mm additional thickness insert may be needed

Alignment rod can be inserted in the gauge to help check alignment.



Gauge



Alignment Rod

PRODUCT NO:
5216-00 [Assembly]
Individual/Replacement Parts:
5216-01 [Gauge]
Overall Length: 7.5" (19,1 cm)
Width: 2.5" (5,4 cm)
Prong Length: 2" (5,1 cm)
Prong Width: (.75" (1,9 cm)
Gap Between Prongs: 1" (2,5 cm)
5216-02 [Alignment Rod]
Overall Length: 18" (45,7 cm)
Diameter: .1875" (4,75 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.