



Taper Head Impactor

Designed by Byron E. Dunaway, MD & Wayne Goldstein, MD

Designed to impact a modular head during minimally invasive THR



The impactor has a protective coating to interface against the implant to help prevent damage while seating the implant. Can be used with 22mm to 32mm heads. Steam and gas sterilizable.

PRODUCT NO:
7840
Overall Length: 12" (30,5 cm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362