

Teurlings Modified Bent Hohmann Retractor

Designed by Luc Teurlings, MD

A twisted hohmann retractor for TKA, can help to protect femoral cuts during MCL/LCL retraction



PRODUCT NO:

7109

Overall Length: 8.625" (21,9 cm)

Depth from Bend: 4.25" (10,8 cm)

Blade Width (Upper Portion): 17 mm

Blade Width (Lower Portion): 11 mm



ISO 13485:2016
FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362