



Tibia AccuAngle

Designed to be placed on the cutting block to check if the cut is level




Product No: 1145

Shown at Actual Size

DIMENSIONS: 3"x 2"



Two magnets in the base help to hold the AccuAngle in place upon a cutting block.



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2007 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000

Fax 912 236 7766

Internet www.innomed.net

eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

