

# Universal Calibrated Tibia/Patella Clamp

Designed by S. David Stulberg, MD



*Designed to be used to remove a tibia wedge, helps in everting the patella, and calibrations help in measuring the thickness of the patella and tibia wedges*

**PRODUCT NO:**

3685

Overall Length: 10" (25,4 cm)  
Calibrations: 0 to 26 mm

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**