



# Uni Medial/Lateral Ligament Retractor

Designed by Kurt Kramer, PA-C

*Designed to be placed in the medial/lateral tibial recess while making the horizontal tibial cut during unicompartmental knee arthroplasty—helping to retract and protect the medial and lateral collateral ligaments*

Ambidextrous, ergonomic design allows for comfortable and natural hand positioning, helping to improve MCL/LCL protection and ease of use, especially in the obese patient.



**PRODUCT NO:**  
**3632**  
 Overall Length: 4.25" (10,8cm)  
 Blade Width: 8.8mm  
 Blade Depth: 2.375" (6cm)



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
 www.innomed.net info@innomed.net

912.236.0000 Phone  
 912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
 Fax +41 41 740 67 71



**1.800.548.2362**