




Vaughan Distal Bicep Tendon Repair Retractor

Designed by Roderick A. Vaughan, MD

Designed to retract in a continuous way in three directions, helping to prevent the surrounding vital structures from entering the field while drilling or performing the repair work

PRODUCT NO:

3223
Overall Length: 8.375" (21,3 cm)
Handle Length: 5.25" (13,3 cm)
Depth: 2" (5,1 cm)



Profile View

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED