


Before

After



# Vaughan Endzone Retractor

Designed by Roderick Vaughan, MD

*Designed for use when placing the end screws while plating a fracture using a minimally invasive technique*

The "U"-shaped wall design helps allow the maximal exposure along the length, or "endzone", of an incision while maintaining adequate width and retraction along the sides of the exposure.



**PRODUCT NO:**

1766

Overall Length: 8.75" (22,2 cm)

Deep Depth: 45 mm

Deep Internal Width: 14 mm

Shallow Depth: 25 mm

Shallow Internal Width: 12 mm



USA MADE

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

# INNOMED

**1.800.548.2362**

© 2017 Innomed, Inc.