



# Vosburg Cannulated Periarticular Clamp

Designed by Caleb Vosburg, MD

*Cannulated clamp tips allow passage of pins up to 7/64" (2.8 mm)*

By compressing the fracture with the clamp and then passing two pins, the clamp can then be removed to allow more working room and versatility when applying a plate.

**PRODUCT NO:**

1864

Overall Length: 13" (33 cm)  
Handle Length: 8" (20.3 cm)  
Ratcheted Opening from 2" to 3.5" (5,4 to 7,6 cm)  
Accepts Pins up to: 7/64" (2.8 mm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**