




# Weinraub Joint and Calcaneal Spreader

Designed by Glenn M. Weinraub DPM, FACFAS

*Designed to assist in the opening of small joints of the foot and hand for the application of fusion and graft techniques*

Provides excellent joint exposure without blocking intra-articular or osteotomy access. Helps prevent slippage or falling out of the joint by placing the arms on either side of the area to be distracted, driving two pins and opening the joint.



**ORIGINAL DESIGN!**

**New!**

Speed lock helps allow precise control and prevent unintended release

PRODUCT NO'S:	
1870 [Standard 1.6 mm]	<i>New!</i>
Overall Length: 7" (17,8 cm) Pin Diameter: Up to .062" (1/16") (1.6 mm)	
1872 [Standard 2.8 mm]	<i>New!</i>
Overall Length: 7" (17,8 cm) Pin Diameter: Up to .11" (7/64") (2.8 mm)	
1870-SL [Speed Lock 1.6 mm]	<i>New!</i>
Overall Length: 7" (17,8 cm) Pin Diameter: Up to .062" (1/16") (1.6 mm)	
1872-SL [Speed Lock 2.8 mm]	<i>New!</i>
Overall Length: 7" (17,8 cm) Pin Diameter: Up to .11" (7/64") (2.8 mm)	

MADE EXCLUSIVELY FOR INNOMED IN GERMANY

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to Launch Our Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

© 2021 Innomed, Inc.