


# Wells Modified Lambotte PAO Osteotomes

Designed by Joel Wells, MD

*Designed to focus on the posterior column osteotomy and connection to the ischial cut – straight, curved and two offset options helps the posterior column osteotomy to be cut with more control*

Silicone handle designed for better control.



| PRODUCT NO'S:                                 |   |
|---|---|
| <b>5276-00</b>                                | [Set of 1 each with case]   |
| <b>Set Includes / Available Individually:</b> |   |
| <b>5276-01</b>                                | [Straight PAO Osteotome]<br>Overall Length: 14" (35,6 cm)<br>Handle Length: 4.5" (11,4 cm)<br>Osteotome Width: 11,1 mm          |
| <b>5276-02</b>                                | [Curved PAO Osteotome]<br>Overall Length: 13.875" (35,3 cm)<br>Handle Length: 4.5" (11,4 cm)<br>Osteotome Width: 12,7 mm        |
| <b>5276-03</b>                                | [Offset PAO Osteotome]<br>Overall Length: 13.625" (34,6 cm)<br>Handle Length: 4.5" (11,4 cm)<br>Osteotome Width: 12,7 mm        |
| <b>5276-04</b>                                | [Offset Curved PAO Osteotome]<br>Overall Length: 13.625" (34,6 cm)<br>Handle Length: 4.5" (11,4 cm)<br>Osteotome Width: 12,7 mm |
| <b>9007</b>                                   | [Case]  |

MADE EXCLUSIVELY  
FOR INNOVATED BY  
GERMANY

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362