

Whang Tibial Osteotome

Designed by William Whang, MD

Designed to disrupt the interface of a well fixed tibial base, specifically the lateral portion



PRODUCT NO:

5338

Overall Length: 8" (20,3 cm)

Handle Length: 4.5" (11,4 cm)

Blade Thickness: 2,5 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362