

# Whelan Femoral Neck Elevator

Designed by E. J. Whelan, III, MD

*Elevator has long tines to rest on the stronger bone at the base of the neck and calcar, and also fits well over the lesser trochanter and iliopsoas tendon for femoral broaching*



**PRODUCT NO:**

3414

Overall Length: 13.75" (34,9 cm)

Depth from Bend: 1.5" (3,8 cm)

Blade Width: 2,4 cm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**