



Wiater Shoulder Slide

Designed by J. Michael Wiater, MD, FAAOS, FAOA

Designed to help avoid damage to the prosthetic bearing surfaces during dislocation and reduction of a shoulder arthroplasty

Also useful for total hip arthroplasty or hip preservation procedures in smaller patients.

Manufactured of delrin to help eliminate damage to the implant. Can be steam or gas sterilized and is radiolucent.

PRODUCT NO:

6879

Overall Length: 11" (27,9 cm)
Width: 1.375" (3,5 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED 