



Sharp and dull teeth available in two depths

# Wide Rake Retractors with Ergonomic Handle

*Designed for general use soft tissue retraction, the ergonomic handle allows for a better grip and less fatigue*

Non-glare finish featured on the metal retractor parts.



PRODUCT NO'S:	
6051 [Deep, Sharp] Overall Length: 11.175 (28,3 cm) Blade Width: 2.375" (6 cm) Blade Depth: 2.75" (7 cm)	6053 [Shallow, Sharp] Overall Length: 11.175 (28,3 cm) Blade Width: 2.375" (6 cm) Blade Depth: 1.875" (4,8 cm)
6052 [Deep, Blunt] Overall Length: 11.175 (28,3 cm) Blade Width: 2.375" (6 cm) Blade Depth: 2.75" (7 cm)	6054 [Shallow, Blunt] Overall Length: 11.175 (28,3 cm) Blade Width: 2.375" (6 cm) Blade Depth: 1.875" (4,8 cm)



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.