

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

**Innomed**  
Instrument Division  
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Woolley Tibia Punch

The Woolley Tibia Punch is designed to impact cancellous bone in the subchondral weight bearing region of the tibia. This helps to improve the mechanical interlock in the cancellous bone/cement interface. The sharp tips can be used on normal and dense cancellous bone, and they can also be used when a significant deformity has been encountered resulting in sclerotic bone.

*Designed by*  
D. Woolley, M.D.

Product No: 5140

*Designed to impact cancellous bone to help improve bone/cement interface*

