




Extra Long Rongeur

*Helpful in minimally invasive total hip surgery
by keeping hands out of the field of view*



Product No's:
Overall Length: 14"
1771-01 5 x 16mm Jaw Bite
1771-02 8 x 16mm Jaw Bite
1771-03 12 x 16mm Jaw Bite



Mazzara Pistol Grip Extra Long Rongeur

Designed by
James T. Mazzara, MD

Product No: 1768-02
Jaw Bite: 8 x 16mm
Overall Length: 12.5"
Shaft-to-End Length: 6"

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

