



Yezerksi Small Bone Rongeur

Designed by John Yezerksi, MD

Designed for small bone applications in the hand and foot



Small

Extra Small

PRODUCT NO'S:

1789 [Small]

Overall Length: 7.125" (18,1 cm)

Jaw Width: 4 mm

Jaw Bite Width: 3 mm

Jaw Bite Length: 20 mm

1789-01 [Extra Small]

Overall Length: 4.5" (11,4 cm)

Jaw Width: Tapers from 4,6 mm to 2 mm

Jaw Bite Length: 11 mm



Scan to
Launch Our
Website



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362