

“Z” Knee Retractor

Helps create better access to the articulating surfaces


Designed to expose the femur and the tibia during knee surgery for better access to the articulating surfaces. The “Z” contouring of the retractor provides the surgeon with an open field of view and working area.

PRODUCT NO:

4420-00

Overall Length: 7.25" (18,4 cm)

Blades: 11 mm Wide, 3" (7,6 cm) Deep



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.