

Zelicof Winged Retractor

Designed by Stephen B. Zelicof, MD Ph.D

A non-modular, fixed point retractor for pelvic visualization in THA and pelvic surgery, with the distal post and wing providing bony stability and retraction




Thin Shaft
#6117-5

Thick Shaft
#6117-10

PRODUCT NO:

6117-5 [Thin Shaft]
 Shaft Diameter: 6 mm
 Overall Length: 9" (22.9 cm)
 Pad Dimensions: 31,4 x 19,6 mm
 Prong Length: .85" (21,6 mm)
 Prong Width (before taper): 3 mm

6117-10 [Thick Shaft]
 Shaft Diameter: 10 mm
 Overall Length: 9" (22,9 cm)
 Pad Dimensions: 31,4 x 19,6 mm
 Prong Length: .85" (21,6 mm)
 Prong Width (before taper): 3 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
 www.innomed.net info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
 Fax +41 41 740 67 71

1.800.548.2362