


# Zell Fixed Angle Wire Guide

Designed by Richard Zell, MD

*Designed to help with placement of guide wires for cannulated screws and k-wires in foot and ankle surgeries, such as bunion surgery, midfoot fusion, and midfoot ORIF*



**PRODUCT NO:**

3021

Overall Length: 2.75" (7 cm)

Handle Platform: 1" x .875" (25 x 22 mm)

Guide Tube Angle: 35°

Guide Tube for wires up to: .052"/1.3 mm



ISO 13485:2016  
**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**