

STERILIZATION PROCEDURE

Stainless Steel Instruments

Please use the following guidelines when sterilizing this product:

Sterilizer Type:	Pre-vacuum Steam Sterilizer
Minimum Temperature:	134 degrees C / 270 degrees F
Minimum Cycle Time:	3 minutes
Minimum Dry Time:	10 minutes

Sterilizer Type:	Flash Sterilizer
Minimum Temperature:	134 degrees C / 270 degrees F
Minimum Cycle Time:	3 minutes
Minimum Dry Time:	0 minutes

Revised 1/1/12

FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

1.800.548.2362