





Knee Retractor System

Designed by S. David Stulberg, MD

Helps free assisting personnel while providing excellent exposure

The Knee Retractor System holds retractors utilizing Velcro® straps. This helps eliminate obstruction of the surgeon's operative area and frees assisting personnel. Four retractor styles are available; straps are available in two lengths. Retractors and straps are autoclavable. The Retractors can be used singularly or in combination.

Velcro® is a registered trademark of Velcro U.S.A.



MIS Modified Wide PCL Retractor with Strap

PRODUCT NO:
3515
Overall Length: 10"
Blade Width Above Prongs: 34mm
Prong Width: 8.5mm | 17mm Gap | 8.5mm



Wide PCL Retractor with Strap

PRODUCT NO:
3525
Overall Length: 10"
Blade Width Above Prongs: 57mm
Prong Width: 8.5mm | 17mm Gap | 8.5mm

Single Prong Collateral Ligament Retractor with Strap

PRODUCT NO:
6650
Overall Length: 8.25"
Blade Width: 14mm



Long Prong Collateral Ligament Retractor with Strap

PRODUCT NO:
6630
Overall Length: 8"
Overall Blade Width: 21mm
Prong Width: 4.5mm | 12mm Gap | 4.5mm



Designed by B. Stubbs, MD

Stubbs Short Prong Collateral Ligament Retractor with Strap

PRODUCT NO:
6640
Overall Length: 8"
Blade Width Above Prongs: 27mm
Prong Width: 4.8mm | 3.4mm Gap | 4.8mm



Strap

PRODUCT NO'S:
Packages of 10
8100-P [Long Strap-Femur]
8120-P [Short Strap-Tibia]

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362