



Amstutz Hip Resurfacing Instruments

Designed by Harlan C. Amstutz, MD

Amstutz Femoral Head-Neck Elevator

Used for acetabular exposure



PRODUCT NO'S:	
3410 [Wide Blade]	Overall Length: 12.25" (31,1 cm) Blade Width at Widest: 67 mm
3410-01 [Narrow Blade]	Overall Length: 12.25" (31,1 cm) Blade Width at Widest: 50 mm



Amstutz Charnley-type Acetabular Exposure Pin Set

Used to enhance exposure in the acetabulum



PRODUCT NO'S:	
1200-00 [Set: Ins/Ext & Two Pins]	
1200-0A [Set: Ins/Ext & Two Pins w/Stop]	
Also sold Individually:	
1200-01 [Inserter/Extractor]	
1200-02 [Pin]	Overall Length: 4.5" (11,4 cm) Pin Depth: 2" (5,1 cm) Pin Diam.: 3.9 mm
1200-03 [Pin with Stop]	Overall Length: 4.5" (11,4 cm) Pin Tip-to-Stop Depth: .75" (1,9 cm) Pin Diameter: 3.2 mm
1200-04 [Deep Pin]	Overall Length: 7" (17,8 cm) Pin Depth: 4.5" (11,4 cm) Pin Diameter: 3.9 mm

Cement Packer & Trimmer



PRODUCT NO:	
4995	Overall Length: 9.75" (24,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

