


# Amstutz Charnley-type Acetabular Exposure Pin Set

Designed by Harlan C. Amstutz, MD

*Used to enhance exposure in the acetabulum*



Deep Pin  
#1200-04

Pin  
#1200-02

Pin with Stop  
#1200-03

Inserter/Extractor  
#1200-01

**PRODUCT NO'S:**

1200-00 [Set: Ins/Ext & Two Pins]

1200-0A [Set: Ins/Ext & Two Pins w/Stop]

**Also sold Individually:**

1200-01 [Inserter/Extractor]

1200-02 [Pin]

Overall Length: 4.5" (11,4 cm)

Pin Depth: 2" (5,1 cm)

Pin Diam.: 3.9 mm

1200-03 [Pin with Stop]

Overall Length: 4.5" (11,4 cm)

Pin Tip-to-Stop Depth: .75" (1,9 cm)

Pin Diameter: 3.2 mm

1200-04 [Deep Pin]

Overall Length: 7" (17,8 cm)

Pin Depth: 4.5" (11,4 cm)

Pin Diameter: 3.9 mm



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**