


Andrews Modified Tibial Wedge Clamp

Designed by Scott Andrews, MD and Kuldeep Sidhu, MD

Designed to help remove the cut tibial bone during total knee procedures

The bone is held by the spikes which helps it to come out in one piece, and also helps with release of soft tissues from the bone.



PRODUCT NO:

3642

Overall Length: 10.25" (26 cm)

Pads: 60 mm x 30 mm

Front Spike Length: 14 mm

Back Spike Length: 7,5 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362