


# Andrews Modified Tibial Fragment Grasper

Designed by Scott Andrews, MD

*Designed to help remove tibial bone during unicondylar and total knee arthroplasty*



**PRODUCT NO:**

**1721**

Overall Length: 10" (25,4 cm)  
Jaw Dimensions: 1.44" x .72" (36,6 x 18,3 mm)  
Lower Jaw Thickness: 1 mm

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**